


Adaptation spéciale pour lit :
Special adaptation for bed:

MMO • Tempo

C-Tempo • 4655003489

MMO®




| C-Tempo | 4655003489 | Montage complet | Complete set |
|------------|------------|-------------------------------------|---|
| | | Le cadre est composé de : | The frame consists of these parts: |
| PRT Tempo | 4655003471 | Portique Tête de lit MMO Tempo (x1) | Head porch bar for MMO Tempo (x1) |
| PRP Tempo | 4655003463 | Portique Pied de lit MMO Tempo (x1) | Foot porch bar for MMO Tempo (x1) |
| B 210 | 4655000055 | Barre octogonale de 210cm (x1) | Plain bar of 210cm (71") (x1) |
| DME Equi* | 4655003059 | Double mâchoire (x2) | Double clamp (x2) |
| T 40 Equi* | 4655003067 | Tube à mâchoire fixe de 40cm (x1) | Fixed clamp bar of 40cm (16") (x1) |
| P 753 | 4655000948 | Etrier de traction (x1) | Handgrip (x1) |
| SP 01 | 4635001304 | Sphère de protection (x3) | Wall protection cushion (x3) |
| PPO | 4655000617 | Porte poulie ouverte (x7) | Open pulley with clamp (x7) |


Construit à partir d'un matériau qui allie la résistance et la légèreté, le Système de Traction KINETEC® est d'un emploi facile pour tout le personnel soignant.

Designed from a material which combines resistance and lightness, the KINETEC® TRACTION System is easy to use for all medical staff.

Kinetec SAS garantit son Système de Traction, 2 ans contre tous vices de fabrication, à partir de la date d'achat par le consommateur.

Kinetec SAS guarantees its Traction System 2 years against all defects of manufacture, starting from the date of purchase by the consumer.

 Dispositif médical conforme aux exigences essentielles de la directive 93/42/CEE relative aux dispositifs médicaux de classe I.

 Medical device conforms to the essential requirements of the directive 93/42/CEE relating to the medical devices of class I.

Demandez notre catalogue complet d'accessoires et de kits de traction (poids, corde, attelle de positionnement, etc...)

Ask for our complete catalogue of accessories and traction kits (weights, cord, positioning splints, and more...)