

# soft traction weights

orthopaedic traction



kinetec®

## traction weights

Kinetec Soft Weights are designed for use with orthopaedic traction.

## eco friendly

Kinetec Soft Weights are made using recycled crumbs from hand therapy thermoplastics, helping to reduce plastic going to landfill.

## safer than hard weights

The weights are soft which protects the nursing staff and visitors from painful knocks if they are accidentally kicked or dropped. They are brightly coloured, helping to prevent accidental collisions which are painful to the patient and staff.

## key information

- Range of bright colours for easy identification
- Easy to transport with black nylon carry handles
- Stainless steel circular attachment ring allows correct positioning of the weight
- Can be combined with other Kinetec Soft Weights to achieve correct traction weight

**use:** Attach traction cord to the steel ring to ensure correct direction of pull

**contents:** Recycled plastic crumbs and sand

**care:** Wipe clean with standard antibacterial wipes

CE Medical Class 1

UK Distributor: Kinetec Medical Products UK Ltd,  
7 Redan Hill Estate, Redan Rd, Aldershot, GU12 4SJ



# order information

Please see order details below for product codes and descriptions.

4655005641 - Kinetec Soft Weight -  
0.25kg (Sky Blue)

4655005659 - Kinetec Soft Weight - 0.5kg  
(Yellow)

4655005607 - Kinetec Soft Weight - 1kg  
(Bright Green)

4655005675 - Kinetec Soft Weight - 2kg  
(Red)

4655005683 - Kinetec Soft Weight - 3kg  
(Hot Pink)

4655005691 - Kinetec Soft Weight - 4kg  
(Purple)

**Kinetec Medical Products UK Ltd**

7 Redan Hill Industrial Estate

Redan Road, Aldershot

GU12 4SJ

Tel: 01252 931980

Fax: 01252 329415

[enquiries@kinetecuk.com](mailto:enquiries@kinetecuk.com)

[www.kinetecuk.com](http://www.kinetecuk.com)